

Flexo-Technic Handels GmbH

USED CODE 12130

Solvent recovery Renzmann M 400



Vessel in mild steel.

For all alternative solvents.

Very easy to operate.

Saves money with the first start.

Keeps your plate-washer clean.

Renzmann M 400 with vacuum
Year of make 1999 – good conditions - checked by Renzmann
With ATEX 95 Certificat
Comes with suction pipe and flexible hose to connect with the filling pipe of the distillation vessel. Comes with emptying hose to connect with the vaccum vessel.

PRICE EXCLUDES:


The costs of delivery, insurance, unloading and positioning, installation, training, commissioning and all taxes.

An estimate will be provided on request, should the customer require assistance with the installation of the machine at the customer's premises

DELIVERY TIME:

Delivery could take place 6 weeks from receipt of your firm order, along with the agreed initial payment and clarification of all technical and commercial details.

TERMS OF PAYMENT: 50% with your order - 50% before crating.



Technical information for the *Renzmann M 400*

Unit for the recovery of high boiling solvent.



Year of manufacture	1999
Volume of distillation vessel	640 litre
Constant filling volume	400 litre
Vacuum	Yes
Inline production	Possible, see schematic
Output per hour	60 - 100 litres*
Operation hours	24 hours possible
Heating up time	1 hour
Type of protection	Exe II T 3
Operation temperature	max. 200
Vacuum abs.	max. 50 mbar
Power requirement	400 V / 50 Hz / 3Ph
Power consumption oil-heater	18 kW
Consumption cooling water	0,8 m ³ /h

* Depends on type of solvent, type and percentage of contamination

Manufactured with an ATEX Quality Assurance System Certification

Technical changes are subject to alteration.

15.05.2008 - M 400 - Code 12130